

Mouse Adiponectin protein (active)

Cat. No. GTX48319-pro**Applications** Functional Assay**Species** Mouse**Package**

25 µg

Applications

Application Note

Biological activity was determined by its ability to inhibit the proliferation of murine myeloid cell lines M1. The ED₅₀ for this effect is < 15 µg/ml.

Properties

Form Lyophilized powder**Buffer** Batch dependent (please contact us for details)**Storage** Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.**Region/Sequence** MKGEPGEAAY MYRSAFSVGL ETRVTPNVP IRFTKIFYNQ QNHYDGSTGK FYCNIPLGY YFSYHITVYMK DVKVSFLFKKD KAVLFTYDQY QEKNVDQASG SVLLHLEVG DQVWLQVYGDG DHNGLYADNV NDSTFTGFL YHDTN**Expression System** E. coli**Purity** > 98% by SDS-PAGE and HPLC.**Endotoxin** < 1 EU/µg**Conjugation** Unconjugated**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).