

Human Adiponectin protein (active)

Cat. No. GTX48328-pro**Applications** Functional Assay**Species** Human**Package**

25 µg

Applications

Application Note

Biological activity was determined by a cytotoxicity assay using M1 cells. The expected ED₅₀ for this effect is 0.5-1.0 µg/ml.

Properties

Form Lyophilized powder**Buffer** Batch dependent (please contact us for details)**Storage** Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.**Region/Sequence** PGAEGPRGFP GIQGRKGEPG EGAYVYRSF SVGLETVVTI PNMPIRFTKI FYNQQNHYDG STGKFHCNIP GLYYFAYHIT
VYMKDVKVSL FKKDKAMLFT YDQYQENNVD QASGSVLLHL EVGDQVWLQV YGEGERNGLY ADNDNDSTFT GFLLYHDTN**Expression System** E. coli**Purity** > 98% by SDS-PAGE and HPLC.**Endotoxin** < 1 EU/µg**Conjugation** Unconjugated**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).