

Human TFF-2 protein (active)

Cat. No. GTX48353-pro

Applications	Functional Assay
Species	Human

Package
20 µg

Applications

Application Note

Assay #1: Biological activity was determined by its ability to chemoattract human monocytes using a concentration range of 1.0-10.0 µg/ml. Assay #2: Biological activity was determined by its ability to chemoattract human MCF-7 cells using a concentration range of 1.0-10.0 ng/ml.

Properties

Form	Lyophilized powder
Buffer	Batch dependent (please contact us for details)
Storage	Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.
Region/Sequence	EKPSPQCQCSR LSPHNRTNCG FPGITSDQCF DNGCCFDSSV TGVPWCFHPL PKQESDQCVM EVSDRRNCGY PGISPEECAS RKCCFSNFIF EVPWCFFPKS VEDCHY
Expression System	E. coli
Purity	> 98% by SDS-PAGE and HPLC.
Endotoxin	< 1 EU/µg
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).