

Human PlGF-1 protein (active)

Cat. No. GTX48384-pro

Application	Functional Assay
Species	Human

Package
25 µg

APPLICATION

Application Note

Biological activity was determined by its ability to chemoattract human monocytes of using a concentration range of 1.0-10.0 ng/ml.

PROPERTIES

Form	Lyophilized powder
Buffer	Batch dependent (please contact us for details)
Storage	Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.
Region/Sequence	MLPAVPPQW ALSAGNGSSE VEVVPFQEVW GRSYCRALER LVDVWSEYPS EVEHMFSPSC VLLRCTGCC GDENLHCVPV ETANVTMQLL KIRSGDRPSY VELTFSQHVR CECRPLREKM KPERCGDAVP RR
Expression System	E. coli
Purity	> 95% by SDS-PAGE and HPLC.
Endotoxin	< 1 EU/µg
Conjugation	Unconjugated
Note	<p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p>



For full product information, images and publications, please visit our [website](#).