

## Human NGF Beta protein (active)

Cat. No. GTX48393-pro

Applications	Functional Assay
Species	Human

Package  
20 µg

## Applications

## Application Note

Biological activity was determined by its ability to stimulate chick E9 DRG neurite outgrowth. The ED<sub>50</sub> is < 1.0 ng/ml, corresponding to a specific activity of > 1 x 10<sup>6</sup> units/mg.

## Properties

Form	Lyophilized powder
Buffer	Batch dependent (please contact us for details)
Storage	Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.
Region/Sequence	SSSHPIFHRG EFSVCDSVSV WVGDKTTATD IKGKEVMVLG EVNINNSVFK QYFFETKCRD PNPVDSGCRG IDSKHWNSYC TTHTFVKAL TMDGKQAAWR FIRIDTACVC VLSRKAVR
Expression System	E. coli
Purity	> 98% by SDS-PAGE and HPLC.
Endotoxin	< 1 EU/µg
Conjugation	Unconjugated

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).