

Human BMP4 protein (active)

Cat. No. GTX48468-pro

Applications Functional Assay

Species Human

Package

5 µg

Applications

Application Note

Biological activity was determined by its ability to induce alkaline phosphatase production by ATDC-5 cells.

Properties

Form Lyophilized powder

Buffer Batch dependent (please contact us for details)

Storage Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

Region/Sequence HHSQRARKKN KNCRRHSLYV DFSDVGWNDW IVAPPGYQAF YCHGDCPFPL ADHLNSTNHA IVQTLVNSVN SSIPKACCVPT
TELSAISMLY LDEYDKVVLK NYQEMVVEGC GCR

Expression System HeLa cells

Purity > 97% by SDS-PAGE and HPLC.

Endotoxin < 1 EU/µg

Conjugation Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).