

Human Apolipoprotein E3 protein

Cat. No. GTX48474-pro
Species Human

Package
500 µg

Applications

Application Note

Properties

Form Lyophilized powder

Buffer Batch dependent (please contact us for details)

Storage Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

Region/Sequence MKVEQAVETE PEPRLRQQTE WQSGQRWELA LGRFWDYLRW VQTLSEQVQE ELLSSQVTQE LRALMDETMK ELKAYKSELE
EQLTPVAEET RARLSKELQA AQARLGADME DVCGRVQYR GEVQAMLGQS TEELRVRLAS HLRKLRKRL RDADDLQKRL
AVYQAGAREG AERGLSAIRE RLGPLVEQGR VRAATVGS LA GQPLQERAQA WGERLRARME EMGSRTRDRL DEVKEQVAEV
RAKLEEQAQQ IRLQAEAFQA RLKSWFEPLV EDMQRQWAGL VEKVQAAVGT SAAPVPSDNH

Expression System E. coli

Purity > 90% by SDS-PAGE and HPLC.

Endotoxin < 1 EU/µg

Conjugation Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


For full product information, images and publications, please visit our [website](#).