

MAT1A + MAT2A antibody (DyLight 488)

Cat. No. GTX48603

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human, Rat, Zebrafish, Bovine, Primate

Package
50 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
--------------------	----------------------

WB	Assay dependent
----	-----------------

Not tested in other applications.

Calculated MW	44 kDa. (Note)
----------------------	----------------------------------

Product Note	We do not recommend use of this product for Mouse samples.
---------------------	------------------------------------------------------------

Properties

Form	Liquid
-------------	--------

Buffer	50mM Sodium Borate
---------------	--------------------

Preservative	0.05% Sodium azide
---------------------	--------------------

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. Protect from light.
----------------	-------------------------------------------------------------------------------------------------------------

Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
----------------------	----------------------------------------------------------------------------------

Immunogen	Synthetic peptide made to an internal portion of the human MAT2A protein (within residues 100-200). [Swiss-Prot# P31153]
------------------	--------------------------------------------------------------------------------------------------------------------------

Purification	Purified by antigen-affinity chromatography
---------------------	---------------------------------------------

Conjugation	DyLight488 Wavelength
--------------------	---------------------------------------

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).