

## MAT1A + MAT2A antibody (DyLight 488)

**Cat. No. GTX48603**

<b>Host</b>	Rabbit
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG
<b>Applications</b>	WB
<b>Reactivity</b>	Human, Rat, Zebrafish, Bovine, Primate

**Package**  
50 µl

## Applications

**Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
--------------------	----------------------

WB	Assay dependent
----	-----------------

Not tested in other applications.

<b>Calculated MW</b>	44 kDa. ( <a href="#">Note</a> )
----------------------	----------------------------------

<b>Product Note</b>	We do not recommend use of this product for Mouse samples.
---------------------	--

## Properties

<b>Form</b>	Liquid
-------------	--------

<b>Buffer</b>	50mM Sodium Borate
---------------	--------------------

<b>Preservative</b>	0.05% Sodium azide
---------------------	--------------------

<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. Protect from light.
----------------	---

<b>Concentration</b>	Batch dependent (Please refer to the vial label for the specific concentration.)
----------------------	--

<b>Immunogen</b>	Synthetic peptide made to an internal portion of the human MAT2A protein (within residues 100-200). [Swiss-Prot# P31153]
------------------	--

<b>Purification</b>	Purified by antigen-affinity chromatography
---------------------	---

<b>Conjugation</b>	DyLight488 <a href="#">Wavelength</a>
--------------------	---------------------------------------

**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).