

MAT1A + MAT2A antibody (DyLight 550)

Cat. No. GTX48605

Host	Rabbit
Clonality	Polyclonal
Isotype	lgG
Applications	WB
Reactivity	Human, Rat, Zebrafish, Bovine, Primate

Package 50 μl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution	
WB	Assay dependent	
Not tested in other applications.		
Calculated MW	44 kDa. (<u>Note</u>)	
Product Note	We do not recommend use of this product for Mouse samples.	

Properties	
Form	Liquid
Buffer	50mM Sodium Borate
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. Protect from light.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide made to an internal portion of the human MAT2A protein (within residues 100-200). [Swiss-Prot# P31153]
Purification	Purified by antigen-affinity chromatography
Conjugation	DyLight550 Wavelength
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 07 Page 1 of 1