

Mlf1 Interacting Protein antibody

Cat. No. GTX48709

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human

Package

50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:2000
IHC-P	Assay dependent
ELISA	1:650000

Not tested in other applications.

Calculated MW 48 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	20mM Potassium Phosphate, 150mM NaCl
Preservative	0.01% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A 418 residue recombinant protein corresponding to the carboxy terminal end of human MLF1IP protein.
Purification	Purified by antigen-affinity chromatography. From serum
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX48709 WB Image

Western blot using GeneTex's affinity purified anti-MLF1IP / PBIP1 antibody shows detection of endogenous MLF1IP protein (a tier of four modified protein bands indicated by the arrowheads) in lysates of Hela cells (- lane). Cells treated with MLF1IP / PBIP1 shRNA (+ lane) show no staining. The identities of the higher and lower molecular weight bands are unknown. Primary antibody was used at 1:1,000.



For full product information, images and publications, please visit our [website](#).