

MIP1 beta antibody

Cat. No. GTX48752

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ELISA
Reactivity	Pig

Package
50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500
ELISA	> 1:5000

Not tested in other applications.

Calculated MW 10 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	20mM Potassium Phosphate, 150mM NaCl
Preservative	0.01% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant protein raised in yeast, corresponding to the 69 amino acids of the mature swine CCL4/MIP1 beta protein.
Purification	Protein A purified From serum
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX48752 WB Image

Western blot using GeneTex's protein-A purified Anti-MIP-1 beta (CCL4) antibody shows detection of recombinant pig MIP-1 beta (CCL4) raised in yeast. The protein was purified and resolved by SDS-PAGE, then transferred to PVDF membrane. Membrane was blocked with 3% BSA , and probed with 4 µg/mL primary antibody overnight at 4°C. After washing, membrane was probed with IRDye800 Conjugated Goat Anti-Rabbit IgG at 1:20,000 for 45 min at room temperature.



GTX48752 WB Image

WB analysis of recombinant pig MIP1 beta protein using GTX48752 MIP1 beta antibody.
Dilution : 4 µg/mL



For full product information, images and publications, please visit our [website](https://www.genetex.com).