

PSCDBP antibody

Cat. No. GTX48859

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IP, ELISA
Reactivity	Human, Mouse, Rat

Package
50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000-1:10000
IP	1:100
ELISA	1:5000-1:20000

Not tested in other applications.

Calculated MW 40 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	20mM Potassium Phosphate, 150mM NaCl
Preservative	0.01% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	85 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide corresponding to a region in the C-terminal portion of mouse Cybr.
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX48859 IP Image

Western blot using GeneTex's anti-Cybr antibody shows detection of recombinant Cybr using anti-HA tag antibody after immunoprecipitation using anti-Cybr. Lane 5 contains lysate from control HEK293 cells. Lane 6 shows IP/WB of HEK293 cells expressing Cybr-HA. Lane 7 and 8 are loaded similarly to lanes 5 and 6 but were immunoprecipitated with anti-HA. Lane 9 contains marker. Lane 10 contains purified CybrHA.



For full product information, images and publications, please visit our [website](#).