

Blocking Buffer, TBS with 1% Casein

Cat. No. GTX48879
Applications

WB, ELISA

Package

250 ml

PRODUCT
Summary

Tris Buffered Saline with 1% casein for blocking of western blots

Applications
Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution
Recommended dilution

WB

Assay dependent

ELISA

Assay dependent

Not tested in other applications.

Properties
Form

Liquid

Storage

Upon receipt store product at 4°C. Product is shipped at ambient temperature.

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


 For full product information, images and publications, please visit our [website](#).