

Blocking Buffer, TBS with 5% Milk

Cat. No. GTX48883

Applications WB, ELISA

References (1)
 Package
 250 ml

PRODUCT

Summary Tris Buffered Saline with 5% low-fat milk for blocking of western blots

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Properties

Form Liquid

Storage Upon receipt store product at 4°C. Product is shipped at ambient temperature.

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
 Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).