

## Blocking Buffer, TBS with 5% Milk

Cat. No. GTX48883

**Application** WB, ELISA**Package**  
250 ml

## PRODUCT

**Summary** Tris Buffered Saline with 5% low-fat milk for blocking of western blots

## APPLICATION

**Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

**Suggested dilution** **Recommended dilution**

WB Assay dependent

ELISA Assay dependent

Not tested in other applications.

## PROPERTIES

**Form** Liquid**Storage** Upon receipt store product at 4°C. Product is shipped at ambient temperature.**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.  
  
Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.For full product information, images and publications, please visit our [website](#).

Date 2024 / 04 / 28 Page 1 of 1