

# ELISA Wash Buffer (30X)

**Cat. No. GTX48885**
**Applications** ELISA

**Package**  
50 ml

## PRODUCT

### Summary

ELISA Wash Buffer is a Tris-based wash buffer for use in cytokine and other sandwich ELISA procedures. The 30X concentrated stock solution is easily diluted 1:30 (vol:vol) in ultrapure water to yield sufficient wash buffer for a typical 96-well microplate ELISA protocol. Supplement an existing ELISA Kit or use as a convenient general purpose wash buffer for customized ELISA methods.

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

### Suggested dilution

### Recommended dilution

ELISA

Assay dependent

Not tested in other applications.

## Properties

### Form

Liquid

### Storage

Store at -20°C immediately upon receipt.

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).