

Trident TBS (tablets)

Cat. No. GTX48886

Applications

WB, ELISA, Cell Culture

Package

100 tablet

PRODUCT

Summary

Pre-measured Tris Buffered Saline Tablets (each sufficient for 100 mL)

Applications

Application Note

When dissolved in 100ml dH₂O, each tablet prepares 100 ml of a 1X TBS solution.

Properties

Form

Liquid

Storage

Store at room temperature

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).