

Trident PBS (tablets)

Cat. No. GTX48887

Application WB, ICC/IF, FACS, ELISA, IHC

Reference (2)

Package

100 tablet

PRODUCT

Summary Pre-measured Phosphate Buffered Saline Tablets (each sufficient for 100 mL)

APPLICATION

Application Note

Each tablet prepares a 1X PBS solution when dissolved in 100 ml of ddH₂O.

PROPERTIES

Storage Store at room temperature.

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).