

## Trident PBS (tablets)

**Cat. No. GTX48887**

**Applications**

WB, ICC/IF, FCM, ELISA, IHC

References ( 2 )

Package

100 tablet

### PRODUCT

**Summary**

Pre-measured Phosphate Buffered Saline Tablets (each sufficient for 100 mL)

### Applications

**Application Note**

Each tablet prepares a 1X PBS solution when dissolved in 100 ml of ddH<sub>2</sub>O.

### Properties

**Storage**

Store at room temperature.

**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).