

PBS (20X)

Cat. No. GTX48888

Applications WB, ICC/IF, FCM, ELISA, IHC

Package
500 ml

PRODUCT

Summary 20X Phosphate Buffered Saline solution

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
--------------------	----------------------

WB	Assay dependent
----	-----------------

ICC/IF	Assay dependent
--------	-----------------

FCM	Assay dependent
-----	-----------------

ELISA	Assay dependent
-------	-----------------

IHC	Assay dependent
-----	-----------------

Not tested in other applications.

Properties

Form Liquid

Storage Store at room temperature.

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).