

# PBS (20X)

**Cat. No. GTX48888**

**Applications** WB, ICC/IF, FCM, ELISA, IHC

**Package**  
500 ml

## PRODUCT

**Summary** 20X Phosphate Buffered Saline solution

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent
FCM	Assay dependent
ELISA	Assay dependent
IHC	Assay dependent

Not tested in other applications.

## Properties

**Form** Liquid

**Storage** Store at room temperature.

**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).