

TBS (20X)

Cat. No. GTX48889

Applications

WB, ICC/IF, FCM, ELISA, IHC

References (1)

Package

500 ml

PRODUCT

Summary

20X Tris Buffered Saline Solution

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution

Recommended dilution

WB Assay dependent

ICC/IF Assay dependent

FCM Assay dependent

ELISA Assay dependent

IHC Assay dependent

Not tested in other applications.

Properties

Form Liquid**Storage** Store at room temperature.

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).