

MYBBP1A antibody, C-term

Cat. No. GTX49078

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Mouse

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.2-2.5 ug/ml
IHC-P	2-10 ug/ml

Not tested in other applications.

Calculated MW 152 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 2% Sucrose
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5-1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide corresponding to a C-terminal region of Mouse MYBBP1A
Purification	Protein A purified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX49078 IHC-P Image

IHC-P analysis of mouse kidney tissue using GTX49078 MYBBP1A antibody at 4.0-8.0µg/ml.


GTX49078 IHC-P Image

IHC-P analysis of mouse lung tissue using GTX49078 MYBBP1A antibody at 4.0-8.0µg/ml.


GTX49078 WB Image

WB analysis of SP2/0 cells using GTX49078 MYBBP1A antibody at 1.25µg/ml.



For full product information, images and publications, please visit our [website](#).