

MEGF11 antibody, Internal

Cat. No. GTX49102

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, IHC
Reactivity	Human, Mouse

References (2)

Package

50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.2-2.5 ug/ml
IHC-P	Assay dependent
IHC	Assay dependent

Not tested in other applications.

Calculated MW 111 kDa. ([Note](#))

Product Note KO/KD validation is based on published data (PMID: 32415115).

Properties

Form	Liquid
Buffer	PBS, 2% Sucrose
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5-1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide corresponding to an Internal region of Human MEGF11
Purification	Affinity Purified
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTx49102 WB Image


WB analysis of human fetal muscle tissue using GTx49102 MEGF11 antibody at 0.2-1 µg/ml.


GTx49102 WB Image

The data was published in the journal Sci Rep in 2020. [PMID: 32415115](https://pubmed.ncbi.nlm.nih.gov/32415115/)


GTx49102 WB Image

The data was published in the journal Sci Rep in 2020. [PMID: 32415115](https://pubmed.ncbi.nlm.nih.gov/32415115/)


GTx49102 WB Image

The data was published in the journal Sci Rep in 2020. [PMID: 32415115](https://pubmed.ncbi.nlm.nih.gov/32415115/)



For full product information, images and publications, please visit our [website](https://www.genetex.com).