

## GTF2H1 antibody [3E8-D3-H8]

Cat. No. GTX49172

|              |            |
|--------------|------------|
| Host         | Mouse      |
| Clonality    | Monoclonal |
| Isotype      | IgG1       |
| Applications | WB, IP     |
| Reactivity   | Human      |

Package  
100 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | Assay dependent      |
| IP                 | Assay dependent      |

Not tested in other applications.

Calculated MW 62 kDa. ( [Note](#) )

## Properties

|               |                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                     |
| Buffer        | PBS, 0.5% BSA, 50% glycerol                                                                                                                                                                                                |
| Preservative  | 0.02% Sodium azide                                                                                                                                                                                                         |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | Batch dependent (Please refer to the vial label for the specific concentration.)                                                                                                                                           |
| Immunogen     | Purified recombinant human TFIIH protein fragments expressed in E.coli.                                                                                                                                                    |
| Purification  | Purified by affinity chromatography                                                                                                                                                                                        |
| Conjugation   | Unconjugated                                                                                                                                                                                                               |


## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

**GTX49172 IP Image**

IP analysis of HeLa lysates using control IgG and GTF2H1 antibody. The precipitates were detected by the same GTF2H1 antibody.


**GTX49172 WB Image**

WB analysis of HeLaNE, Jurkat and HeLa lysates using GTF2H1 antibody at a dilution of 1:300.



For full product information, images and publications, please visit our [website](https://www.genetex.com).