

FXYD2 antibody

Cat. No. GTX49356

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

WB: 1:500 - 1:1000


Calculated MW	7 kDa. (Note)
Product Note	FXYD2 pAb detects endogenous levels of FXYD2 protein.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide corresponding to residues in Human FXYD2.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTG49356 WB Image

Western blot (WB) analysis of FXVD2 pAb at 1:500 dilution Lane1:HepG2 whole cell lysate Lane2:Mouse lung tissue lysate Lane3:Rat brain tissue lysate



For full product information, images and publications, please visit our [website](#).