

REG4 antibody

Cat. No. GTX49362

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

WB: 1:500 - 1:1000

Calculated MW	18 kDa. (Note)
Product Note	REG4 pAb detects endogenous levels of REG4 protein.


Properties

Form	Liquid
Buffer	PBS
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide corresponding to residues in human REG4
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX49362 WB Image

Western blot (WB) analysis of REG4 pAb at 1:500 dilution. Lane1:HEK293T whole cell lysate. Lane2:Raw264.7 whole cell lysate. Lane3:H9C2 whole cell lysate.



For full product information, images and publications, please visit our [website](https://www.genetex.com).