

## Bisbenzimidazole / Hoechst 33342

Cat. No. GTX49386

**Application** ICC/IF, FACS**Package**  
100 mg

## PRODUCT

**Summary**

This product is soluble in water (50 mg/ml), yielding a clear solution. The pH of a 2% solution is 1.9. It has been observed that this material will precipitate from phosphate buffer solutions.

## APPLICATION

**Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

**Suggested dilution****Recommended dilution**

ICC/IF

Assay dependent

FACS

Assay dependent

Not tested in other applications.

## PROPERTIES

**Form**

Lyophilized powder

**Storage**

Aqueous solutions are stable for 1 month if kept in the dark at 2-8°C. Store at -20°C.

**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).