

# Bisbenzimidazole / Hoechst 33342

**Cat. No. GTX49386**

**Applications** ICC/IF, FCM

**Package**  
100 mg

## PRODUCT

### Summary

This product is soluble in water (50 mg/ml), yielding a clear solution. The pH of a 2% solution is 1.9. It has been observed that this material will precipitate from phosphate buffer solutions.

### Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

### Suggested dilution

### Recommended dilution

ICC/IF

Assay dependent

FCM

Assay dependent

Not tested in other applications.

## Properties

### Form

Lyophilized powder

### Storage

Aqueous solutions are stable for 1 month if kept in the dark at 2-8°C. Store at -20°C.

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).