

Bisbenzimidazole / Hoechst 33342

Cat. No. GTX49386
Applications ICC/IF, FCM

Package
100 mg

PRODUCT

Summary

This product is soluble in water (50 mg/ml), yielding a clear solution. The pH of a 2% solution is 1.9. It has been observed that this material will precipitate from phosphate buffer solutions.

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution

Recommended dilution

ICC/IF

Assay dependent

FCM

Assay dependent

Not tested in other applications.

Properties

Form

Lyophilized powder

Storage

Aqueous solutions are stable for 1 month if kept in the dark at 2-8°C. Store at -20°C.

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).