

KSHV ORF73 antibody [AT4C11]

Cat. No. GTX50023

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, IHC-P, ELISA
Reactivity	Kaposi's sarcoma-associated herpesvirus

Package
100 µl

Applications

Application Note

We recommend the following starting dilutions:

Western blot is 1:1,000 ~ 1:2,000


Immunohistochemistry: Use at 1:50~100.

Optimal working concentrations should be determined experimentally by the end user.


Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant human ORF73/HHV8 (122-329aa) purified from E. coli
Purification	By protein-G affinity chromatography
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX50023 WB Image

Transfected and non-transfected 293T cell lysates were resolved by SDS-PAGE, transferred to NC membrane and probed with anti-KSHV ORF73 (1:1,000) and anti-MBP (1:1,000). Proteins were visualized using a goat anti-mouse secondary antibody conjugated to HR


GTX50023 IHC-P Image

Paraffin embedded sections of Kaposi's sarcoma tissue were incubated with anti-KSHV ORF73 (1:50) for 2 hours at room temperature. Antigen retrieval was performed in 0.1M sodium citrate buffer and detected using Diaminobenzidine (DAB)



For full product information, images and publications, please visit our [website](https://www.genetex.com).