

KIR2D antibody [2H6]

Cat. No. GTX50072

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB, ELISA
Reactivity	Human

Package
100 µl

Applications

Application Note

We recommend the following starting dilutions:

Western Blot: Use at 1:500 ~ 1,000.

Optimal working concentrations should be determined experimentally by the end user.


Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant human kIR2D (44-202aa) purified from E. coli
Purification	By protein-G affinity chromatography
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX50072 WB Image

Human recombinant protein kIR2DL1, kIR2DL3 and kIR2DL4 (each 20ng per well) were resolved by SDS-PAGE, transferred to PVDF membrane and probed with anti-human kIR2D (1:500). Proteins were visualized using a goat anti-mouse secondary antibody conjugated t



For full product information, images and publications, please visit our [website](#).