

alpha A Crystallin antibody [c9F2]

Cat. No. GTX50077

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

We recommend the following starting dilutions:

Western blot analysis: 1:500 ~ 1:2000.


Optimal working concentrations should be determined experimentally by the end user.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant human alpha A Crystallin (1-173aa) purified from E. coli
Purification	By protein-G affinity chromatography
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](https://www.genetex.com).

DATA IMAGES

GTX50077 WB Image

Mouse eye extracts and recombinant proteins (Crystallin-alpha A and B) were resolved by electrophoresis, transferred to PVDF membrane and probed with anti-Crystallin alpha A (1:1000). Proteins were visualized using a goat anti-mouse secondary antibody con



For full product information, images and publications, please visit our [website](#).