

KIR2DL3 antibody [190IIC311]

Cat. No. GTX50078

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	ICC/IF, FCM, ELISA
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	Assay dependent
FCM	5-10 µg for 10 ⁶ cells
ELISA	1:1000

Not tested in other applications.

Product Note This antibody is specific for human KIR2DL3 protein, and it does not cross react with human KIR2DL1 and KIR2DL4 protein

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant human KIR2DL3 (19-161aa) purified from E. coli
Purification	By protein-G affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX50078 ICC/IF Image

ICC/IF analysis of HeLa cells using GTX50078 KIR2DL3 antibody [190IIC311].

Green : Primary antibody

Blue : DAPI

Dilution : 1:100



GTX50078 Image

WB analysis of various human recombinant proteins using GTX50078 KIR2DL3 antibody [190IIC311].

Dilution : 1:1000



For full product information, images and publications, please visit our [website](#).