

Estrogen Receptor alpha (phospho Ser167) antibody

Cat. No. GTX50140

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, FCM
Reactivity	Human

References (1)

★★★★☆ Review (1)

Package

100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:1000
ICC/IF	1:100-1:200
IHC-P	1:50-1:100
FCM	Assay dependent

Not tested in other applications.

Calculated MW 66 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 150mM NaCl, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Peptide sequence around phosphorylation site of serine 167 (L-A-S(p)-T-N) derived from human Estrogen Receptor alpha.
Purification	Purified by antigen-affinity chromatography. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Conjugation	Unconjugated

For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 09 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX50140 ICC/IF Image


ICC/IF analysis of methanol-fixed MCF cells using GTX50140 Estrogen Receptor alpha (phospho Ser167) antibody.



GTX50140 WB Image

WB analysis of extracts from MCF-7 cells using GTX50140 Estrogen Receptor alpha (phospho Ser167) antibody.

Lane 1 : Primary antibody pre-incubated with the antigen specific peptide



GTX50140 IHC-P Image

IHC-P analysis of human breast carcinoma tissue using GTX50140 Estrogen Receptor alpha (phospho Ser167) antibody.

Left : Primary antibody

Right : Primary antibody pre-incubated with the antigen specific peptide



For full product information, images and publications, please visit our [website](https://www.genetex.com).