

FOXO1 (phospho Ser256) antibody

Cat. No. GTX50179

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, IP
Reactivity	Human, Mouse

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:1000
ICC/IF	1:100-1:200
IHC-P	1:50-1:100
IP	Assay dependent

Not tested in other applications.

Calculated MW 70 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 150mM NaCl, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Peptide sequence around phosphorylation site of serine 256 (A-A-S(p)-M-D) derived from human FOXO1.
Purification	Purified by antigen-affinity chromatography. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX50179 ICC/IF Image


ICC/IF analysis of methanol-fixed MCF-7 cells using GTX50179 FOXO1 (phospho Ser256) antibody.


GTX50179 WB Image

WB analysis of extracts from 293 cells untreated or treated with serum using GTX50179 FOXO1 (phospho Ser256) antibody.


GTX50179 ICC/IF Image

ICC/IF analysis of methanol-fixed MEF cells untreated or treated with LY2904 using GTX50179 FOXO1 (phospho Ser256) antibody.


GTX50179 IHC-P Image

IHC-P analysis of human breast carcinoma tissue using GTX50179 FOXO1 (phospho Ser256) antibody.

Left : Primary antibody

Right : Primary antibody pre-incubated with the antigen specific peptide



For full product information, images and publications, please visit our [website](#).