

EGFR (phospho Tyr845) antibody

Cat. No. GTX50264

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
--------------------	----------------------

WB	1:500-1:1000
----	--------------

Not tested in other applications.

Calculated MW	134 kDa. (Note)
---------------	-----------------------------------

Product Note

The full-length EGFR includes the 24 aa signal peptide that is cleaved during EGFR maturation. The Tyr845 of mature EGFR is equivalent to Tyr869 of uncleaved EGFR.

Properties

Form	Liquid
------	--------

Buffer	PBS, 150mM NaCl, 50% Glycerol
--------	-------------------------------

Preservative	0.02% Sodium azide
--------------	--------------------

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
---------	--

Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
---------------	--

Immunogen	Peptide sequence around phosphorylation site of tyrosine 845 (K-E-Y(p)-H-A) derived from human EGFR.
-----------	--

Purification	Purified by antigen-affinity chromatography. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
--------------	---

Conjugation	Unconjugated
-------------	--------------

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 03 / 31 Page 1 of 2

DATA IMAGES

**GTX50264 WB Image**

WB analysis of extracts from MDA cells untreated or treated with EGF using GTX50264 EGFR (phospho Tyr845) antibody.



For full product information, images and publications, please visit our [website](#).