

## p62Dok antibody

Cat. No. GTX50588

|              |                   |
|--------------|-------------------|
| Host         | Rabbit            |
| Clonality    | Polyclonal        |
| Isotype      | IgG               |
| Applications | WB, ICC/IF, IHC-P |
| Reactivity   | Human             |

Package  
100 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1:500-1:1000         |
| ICC/IF             | 1:100-1:200          |
| IHC-P              | 1:50-1:100           |

Not tested in other applications.

Calculated MW 52 kDa. ( [Note](#) )

## Properties

|               |                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                     |
| Buffer        | PBS, 150mM NaCl, 50% Glycerol                                                                                                                                                                                              |
| Preservative  | 0.02% Sodium azide                                                                                                                                                                                                         |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | 1 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                   |
| Immunogen     | Peptide sequence around aa.396~400 (P-I-Y-D-E) derived from human p62Dok.                                                                                                                                                  |
| Purification  | Purified by antigen-affinity chromatography.                                                                                                                                                                               |
| Conjugation   | Unconjugated                                                                                                                                                                                                               |

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX50588 WB Image

WB analysis of extracts from K562, HL-60, and JK cells using GTX50588 p62Dok antibody.




### GTX50588 IHC-P Image

IHC-P analysis of human breast carcinoma tissue using GTX50588 p62Dok antibody.

Left : Primary antibody

Right : Primary antibody pre-incubated with the antigen specific peptide



### GTX50588 ICC/IF Image

ICC/IF analysis of methanol-fixed HeLa cells using GTX50588 p62Dok antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).