

## STAT3 antibody

Cat. No. GTX50709

|              |                   |
|--------------|-------------------|
| Host         | Rabbit            |
| Clonality    | Polyclonal        |
| Isotype      | IgG               |
| Applications | WB, ICC/IF, IHC-P |
| Reactivity   | Human             |

Package  
100 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1:500-1:1000         |
| ICC/IF             | 1:100-1:200          |
| IHC-P              | 1:50-1:100           |

Not tested in other applications.

Calculated MW 88 kDa. ( [Note](#) )

## Properties

|               |                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                     |
| Buffer        | PBS, 150mM NaCl, 50% Glycerol                                                                                                                                                                                              |
| Preservative  | 0.02% Sodium azide                                                                                                                                                                                                         |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | 1 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                   |
| Immunogen     | Peptide sequence around aa.703~707(A-P-Y-L-K) derived from human STAT3.                                                                                                                                                    |
| Purification  | Purified by antigen-affinity chromatography.                                                                                                                                                                               |
| Conjugation   | Unconjugated                                                                                                                                                                                                               |

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

## DATA IMAGES




### GTx50709 ICC/IF Image

ICC/IF analysis of methanol-fixed HeLa cells using GTx50709 STAT3 antibody.



### GTx50709 WB Image

WB analysis of extracts from HeLa cells using GTx50709 STAT3 antibody.



### GTx50709 IHC-P Image

IHC-P analysis of human breast carcinoma tissue using GTx50709 STAT3 antibody.



For full product information, images and publications, please visit our [website](#).