

Neurologin 1 + Neurologin 2 antibody

Cat. No. GTX51022

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	FCM, ELISA
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FCM	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Properties

Form	Liquid
Buffer	1% BSA, 50% Glycerol
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide derived from human Neurologin 1 + Neurologin 2.
Purification	Protein A purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX51022 FCM Image

FACS analysis of U251 cells using GTX51022 Neuroligin 1 + Neuroligin 2 antibody.


Green : Primary antibody

Blue : Unstained cells

Orange : Isotype control

Light blue : Secondary antibody only

Dilution : 1:100



GTX51022 ELISA Image

ELISA analysis of 0.2 μ g GTX51022 antigen in 100 μ l buffer using serial diluted GTX51022 Neuroligin 1 + Neuroligin 2 antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).