

ENPP3 antibody

Cat. No. GTX51518

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	IHC-P, FCM
Reactivity	Human, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1:50-400
FCM	1:20-100

Note : Permeabilization : 20% PBST for 20 min at room temperature

Not tested in other applications.

Properties

Form	Liquid
Buffer	1% BSA, 50% Glycerol
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide derived from mouse ENPP3(31-80).
Purification	Protein A purified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTX51518 FCM Image

FACS analysis of HeLa cells using GTX51518 ENPP3 antibody.

Green : Primary antibody

Blue : Isotype control

Dilution : 1:10


GTX51518 FCM Image

FACS analysis of Raji (Positive) and Molt-4 (Negative control) cells using GTX51518 ENPP3 antibody.

Green : Primary antibody

Black : Unstained cells


Orange : Isotype control

Blue : Secondary antibody only

Dilution : 1:50

Fixation : 4% PFA for 10min at room temperature

Permeabilization : 20% PBST for 20 min at room temperature


GTX51518 IHC-P Image

IHC-P analysis of rat brain tissue using GTX51518 ENPP3 antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).