

TACI antibody

Cat. No. GTX51554

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:300-1000
IHC-P	1:50-400

Not tested in other applications.

Calculated MW 27 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	1% BSA, 50% Glycerol
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide derived from mouse TACI(20-70).
Purification	Protein A purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES



GTX51554 IHC-P Image

IHC-P analysis of rat kidney tissue using GTX51554 TACI antibody.

Dilution : 1:200




GTX51554 WB Image

WB analysis of various samples using GTX51554 TACI antibody.

Dilution : 1:200

Lane 1 : mouse brain lysates

Lane 2 : mouse heart lysates



GTX51554 WB Image

WB analysis of mouse spleen tissue lysate using GTX51554 TACI antibody.

Dilution : 1:300



For full product information, images and publications, please visit our [website](https://www.genetex.com).