

SPON1 antibody

Cat. No. GTX51796

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	IHC-P, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1:50-400
ELISA	Assay dependent

Not tested in other applications.

Properties

Form	Liquid
Buffer	1% BSA, 50% Glycerol
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	KLH conjugated synthetic peptide derived from human SPON1(775-807).
Purification	Protein A purified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX51796 IHC-P Image


IHC-P analysis of rat lung tissue using GTX51796 SPON1 antibody.
Dilution : 1:200


GTX51796 IHC-P Image

IHC-P analysis of rat brain tissue using GTX51796 SPON1 antibody.
Antigen retrieval : Boiling in sodium citrate buffer (pH6) for 15min
Dilution : 1:400


GTX51796 IHC-P Image

IHC-P analysis of mouse brain tissue using GTX51796 SPON1 antibody.
Antigen retrieval : Boiling in sodium citrate buffer (pH6) for 15min
Dilution : 1:400


GTX51796 ELISA Image

ELISA analysis of 2µg GTX51796 antigen in 100µl buffer using serial diluted GTX51796 SPON1 antibody.



For full product information, images and publications, please visit our [website](http://www.genetex.com).