

IKK alpha (phospho Ser176) + IKK beta (phospho Ser177) antibody

Cat. No. GTX52310

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC
Reactivity	Human, Mouse, Rat

[References \(3 \)](#)[Package](#)

100 µg

Applications

Application Note

WB 1:500-1:2000, IHC 1:50-1:200

Calculated MW 85 kDa. ([Note](#))**Product Note** CHUK

Properties

Form Liquid**Buffer** PBS, 150mM NaCl, 50% Glycerol**Preservative** 0.02% Sodium azide**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 1 mg/ml (Please refer to the vial label for the specific concentration.)**Immunogen** A synthesized peptide derived from human IKK- beta around the phosphorylation site of Ser177**Purification** Antigen affinity purification**Conjugation** Unconjugated**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).