

ERK1 (phospho Thr202) + ERK2 (phospho Thr185) antibody

Cat. No. GTX52319

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

WB 1:500-1:2000

Calculated MW 43 kDa. ([Note](#))

Product Note MAPK3

Properties

Form Liquid

Buffer PBS, 150mM NaCl, 50% Glycerol

Preservative 0.02% Sodium azide

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 1 mg/ml (Please refer to the vial label for the specific concentration.)

Immunogen A synthesized peptide derived from human p44/42 MAP Kinase around the phosphorylation site of Thr202

Purification Antigen affinity purification

Conjugation Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).