

## AHR (phospho Ser36) antibody

**Cat. No. GTX52324**

|                     |                   |
|---------------------|-------------------|
| <b>Host</b>         | Rabbit            |
| <b>Clonality</b>    | Polyclonal        |
| <b>Isotype</b>      | IgG               |
| <b>Applications</b> | WB, IHC           |
| <b>Reactivity</b>   | Human, Mouse, Rat |

**Package**  
100 µg

## Applications

**Application Note**

WB 1:500-1:2000, IHC 1:50-1:200

**Calculated MW** 96 kDa. ( [Note](#) )**Product Note** AHR

## Properties

|                      |                                                                                                                                                                                                                                                                                                                                             |
|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                                                                                                                                                                                                                                                      |
| <b>Buffer</b>        | PBS, 150mM NaCl, 50% Glycerol                                                                                                                                                                                                                                                                                                               |
| <b>Preservative</b>  | 0.02% Sodium azide                                                                                                                                                                                                                                                                                                                          |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.                                                                                                                  |
| <b>Concentration</b> | 1 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                                                                                                                                    |
| <b>Immunogen</b>     | A synthesized peptide derived from human AhR around the phosphorylation site of Ser36                                                                                                                                                                                                                                                       |
| <b>Purification</b>  | Antigen affinity purification                                                                                                                                                                                                                                                                                                               |
| <b>Conjugation</b>   | Unconjugated                                                                                                                                                                                                                                                                                                                                |
| <b>Note</b>          | <p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p> |



For full product information, images and publications, please visit our [website](#).