

## IL2 Receptor alpha antibody [MAB0817]

**Cat. No. GTX52405**

<b>Host</b>	Rat
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG1
<b>Applications</b>	Depletion, Neutralizing /Inhibition
<b>Reactivity</b>	Mouse

**Package**  
200 µg

## Applications

**Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
Depletion	Assay dependent
Neutralizing /Inhibition	Assay dependent

**Note : Depletion of CD25+CD4+ Treg cells**  
**Growth inhibition of IL-2-dependent T-cell lines**  
**Blocking of IL-2 binding to the low- and high-affinity receptors**

Not tested in other applications.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS
<b>Preservative</b>	No preservatives
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
<b>Concentration</b>	Batch dependent (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	IL-2-dependent cytolytic mouse T-cell
<b>Purification</b>	Protein A/G purified From tissue culture supernatant
<b>Conjugation</b>	Unconjugated

**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).