

## CD28 antibody [MAB0818]

Cat. No. GTX52406

Host	Golden Syrian Hamster
Clonality	Monoclonal
Isotype	IgG
Applications	Activation, Neutralizing/Inhibition
Reactivity	Mouse

Package  
200 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
Activation	Assay dependent
Neutralizing/Inhibition	Assay dependent

**Note : Induction of thymocyte differentiation****Blocking of allogeneic mixed leukocyte response and inhibition of MHC-unrestricted CTL cytotoxicity**

Not tested in other applications.

## Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Mouse T-cell lymphoma
Purification	Protein A/G purified From tissue culture supernatant
Conjugation	Unconjugated

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).