

CD45R antibody [MAB0822]

Cat. No. GTX52410

Host	Rat
Clonality	Monoclonal
Isotype	lgG2
Applications	WB, FCM, IP, Neutralizing/Inhibition
Reactivity	Mouse

Package $200\ \mu g$

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
FCM	Assay dependent
IP	Assay dependent
Neutralizing/Inhibition	Assay dependent
Note: Blocking T cell activation in vivo	

Not tested in other applications.

Product Note This antibody was selected for its ability to detect human CD45RB on most leukocytes.

Properties	
Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Activated mouse B lymphoblasts
Purification	Protein A/G purified From tissue culture supernatant
Conjugation	Unconjugated



For full product information, images and publications, please visit our website.

Date 2025 / 10 / 25 Page 1 of 2



For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 10 / 25 Page 2 of 2